USER MANUAL

Bedienungsanleitung
Instrukcja obsługi
Návod k použití
Manuel d'utilisation
Istruzioni per l'uso
Manual de instrucciones
Használati útmutató
Brugsanvisning
Käyttöohje
Gebruiksaanwijzing
Bruksanvisning
Instruções de utilização
Používateľská príručka
Ръководство за употреба
Οδηγίες χρήσης
Upute za uporabu
Naudojimo instrukcija
Manual de utilizare
Navodila za uporabo

MICROSCOPE

DE	Produktname	Mikroskop
EN	Product name	Microscope
PL	Nazwa produktu	Mikroskop
CZ	Název výrobku	Mikroskop
FR	Nom du produit	Microscope
IT	Nome del prodotto	Microscopio
ES	Nombre del producto	Microscopio
HU	Termék neve	Mikroszkóp
DA	Produktnavn	Mikroskop
FI	Tuotteen nimi	Mikroskooppi
NL	Productnaam	Microscoop
NO	Produktnavn	Mikroskop
SE	Produktnamn	Mikroskop
PT	Nome do produto	Microscópio
SK	Názov produktu	Mikroskop
BG	Име на продукта	Микроскоп
EL	Όνομα προϊόντος	Μικροσκόπιο
HR	Naziv proizvoda	Mikroskop
LT	Produkto pavadinimas	Mikroskopas
RO	Numele produsului	Microscop
SL	Ime izdelka	mikroskop
DE Modell EN Product model PL Model produktu CZ Model výrobku FR Modèle IT Modello ES Modelo HU Modell DA Model FI Tuotteen malli NL Productmodel NO Produktmodell SE Produktmodell PT Modelo do produto SK Model BG Модел на продукт EL Μοντέλο προϊόντος HR Model proizvoda LT : Gaminio modelis RO : Model de produs SL : Model izdelka		SBS-MK-5 SBS-MK-7
DE Hersteller EN Manufacturer PL Producent CZ Výrobce FR Fabricant IT Produttore ES Fabricante HU Termelő DA Producent FI Valmistaja NL Producent NO Produsent SE Tillverkare PT Fabricante SK Výrobca BG Производител EL Κατασκευαστής HR Proizvođač LT Gamintojas RO Producător SL Proizvajalec		expondo Polska sp. z o.o. sp. k.
DE Anschrift des Herstellers EN Manufacturer Address PL Adres producenta CZ Adresa výrobce FR Adresse du fabricant IT Indirizzo del produttore ES Dirección del fabricante HU A gyártó címe DA Producentens adresse FI Valmistajan osoite NL Adres producent NO Produsentens adresse SE Tillverkarens adress PT Endereço do fabricante SK Adresa výrobci BG Адрес на производителя EL : Διεύθυνση κατασκευαστή HR Adresa proizvođača LT Gamintojo adresas RO Adresa producătorului SL Naslov proizvajalca		ul. Nowy Kisielin – Innowacyjna 7, 66-002 Zielona Góra Poland, EU

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- 10- Giaphram
- 11- Fokussierknopf
- 12- Moverule
- 13- Hauptteil

This User Manual has been translated using machine translation. We have made every effort to ensure the translation is accurate, but please note that automated translations are not perfect and are not meant to replace human translators. The official version of the User Manual is in English. Any differences between the translated version and the original English are not legally binding. If you have any questions about the accuracy of the translation, please refer to the English version, which is the official reference. More language versions are available upon request via info@expondo.com.

Technical data

Parameter description	Parameter value	
Product name	Microscope	
Model	SBS-MK-5	SBS-MK-7
Voltage [V]	5	
Magnification ratio	40X	
Dimensions (Width x Length x Height) (mm)	100 x 150 x 190	135 x 180 x 290
Weight [kg]	1.05	1.15
Protection class	II	

Purpose

The product is used to magnify small objects or details that are too small to be seen clearly with the naked eye. By enlarging these objects, the product allows users to observe fine structures and details, making it an essential tool in scientific research, biology, medicine, material science, and various other fields.

Product overview

Assembly

1. Remove the product from the carton and place the foam container on its side.

2. Carefully remove the tape around the perimeter of the foam container and lift the top half.
3. Open the container carefully to prevent the optical items (objectives and eyepieces) from falling and getting damaged.
4. Using both hands (one hand around the arm and the other supporting the base), lift the product out of the container and place it on a stable desk.
5. Take the objectives out of their individual protective vials and install them into the product's nosepiece in order of increasing magnification, turning clockwise from the rear.
6. Insert the eyepiece into the eye tube. The eyepiece is included with the product.

Usage

1. Adjusting the Observation Head:

The observation head rotates 360°, allowing easy adjustment to a comfortable viewing position.

2. Placing the Specimen:

Place the specimen on a glass slide and secure it with the clips or slide-holders of the movable specimen-holder. Ensure the specimen is centered over the stage opening.

3. Illuminating the Specimen:

- For products with a mirror, adjust the mirror to light the entire field of view.
- For products with a built-in electrical illuminator, plug the power cable into the socket and turn on the power switch to illuminate the specimen.

NOTE

Changing the Lamp: Before replacing the lamp, unplug the product and wait for it to cool down to avoid burns. Unscrew the lamp holder, then replace the lamp housing on the base. Mirrors can be changed in the same manner.

4. Interpupillary and Diopter Adjustment (for Binocular Microscopes):

Firmly hold each barrel of the binocular head and adjust the distance between them until you see one circle of light. After interpupillary adjustment, check the graduation value in the middle of the headcover. Rotate the diopter rings to align the same graduation value with the lines on the eye tubes.

5. Focusing the Specimen:

- Rotate the nosepiece to position the lowest power objective over the specimen. A distinct "click" indicates proper alignment.
- Start with the lowest power objective to observe general structural details with a large field of view, then increase magnification as needed for finer details.
- To calculate magnification, multiply the eyepiece power by the objective power.

- Use the coarse focusing knob to move the stage up until the specimen is in focus, then lock the focus-stop screw.
- Be careful to avoid contact between the objective and the slide to prevent damage.
- Use the fine focusing knob to sharpen the image. Look through the eyepiece with either eye, keeping both eyes open. For binocular microscopes, adjust the diopter rings for a sharp, clear image in both eyes.

6. Adjusting the Diaphragm:

- For products with a disc diaphragm, rotate the diaphragm to select an appropriate aperture for suitable background brightness.
- For products with an iris diaphragm, adjust the aperture to control the background brightness.

Maintenance

Always:

- Maintain the product in an environment with the following conditions:
 - **Indoor temperature:** 0°C to 40°C
 - **Maximum relative humidity:** 85%
- Keep the product free from dust and shocks during use.
- Turn off the light immediately after use.
- Clean the optics with a soft lens tissue.
- Use a gentle cleaning solution to clean other parts of the product.
- After use, cover the product with the provided dust cover and store it in a dry, clean place.

Never:


- Wipe the surface of any optics with your fingers.
- Disassemble the objectives or eyepieces for cleaning.
- Mishandle or apply unnecessary force to the product.
- Clean the product with volatile solvents or abrasive cleaners.
- Attempt to service or repair the product yourself.

Disposing of Used Devices

Do not dispose of this device in municipal waste systems. Hand it over to an electric and electrical device recycling and collection point. Check the symbol on the product, instruction manual and packaging. The plastics used to construct the device can be recycled in accordance with their markings. By choosing to recycle you are making a significant contribution to the protection of our environment.

Contact local authorities for information on your local recycling facility.

Parts

- 1- Eyepiece
- 2- Eye tube
- 3- Monocular head
- 4- nosepiece
- 5- objective
- 6- stage
- 7- LED
- 8- Base
- 9- Power socket
- 10- Giaphram
- 11- Focusing-knob
- 12- Moverule
- 13- Main body

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